


# cobas® pro integrated solutions

## Site Requirements V2s1

Sample Supply unit	ISE Unit	cobas c 503 unit	cobas e 801 unit
Up to <b>1000</b> samples per hour processing capacity Up to <b>300</b> samples direct loading Up to <b>300</b> samples direct unloading or bi-directional connectivity to lab automation	Up to <b>300</b> samples per hour Up to <b>900</b> tests per hour	Up to <b>1000</b> tests per hour Reaction time : 3-10 minutes Up to <b>60</b> reagent positions RFID label	Up to 300 tests per hour Up to 48 reagent positions Reaction time: 9/18/27

### Disclaimer



This Pre-Installation Site Checklist does not replace any user or service manual.

The site requirement document can be used as a basis for the pre-installation discussion with the customer.

### Dimensions:

	SSU Including monitor	ISE module	c 503 including sample buffer	e801 Including sample buffer
□ W x D x H	52 x 38 x 56 in.	18 x 46 x 46 in.	61 x 46 x 53 in	61 x 46 x 53 in
□ Weight	551 lbs.	337 lbs.	1336 lbs.	1448 lbs

□ The floor should have a grade (rise or fall) less than 1/200 within the footprint of the system.

### Electrical Power Requirements

- Source
- Two 115 VAC +/- 10% / 60Hz -standard wall receptacle [5-15R 50 AMP (Load tested) - Two (2) options: **United States**
  - connection]
    1. Wall mounted Hubbell # 3771, 3-Wire 50 A, 250VDC - 600VAC
    2. Cord mounted Hubbell# 3762C, 3 Wire 50 A, 250VDC - 600VAC
  - 208 VAC 50 amp 1 connection required (**2 modules**)
  - 208 VAC 50 amp 2 connections required (**3-4 modules**)

Two 15 Amp circuits will supply power to:

1. Roche Service use.
2. HP LaserJet Printer



**Note:** All of the above are Locking Single Receptacles


50 AMP (Load tested) – Two (2) options: **Canada**

1. Wall mounted Hubbell #6369CR, Outlet twist lock non-NEMA 50A
2. Cord mounted Hubbell # CS6364, 50A, 125/250VAC Female Power Plt

**Note:** All of the above are Locking Single Receptacles

### Electrical

One (single phase) 208V 50 Amp receptacle  
Must be located close to the 10 kVA UPS

US & Puerto Rico	Canada
 Wall Receptacle <b>Hubbell 3771</b> (Or Inline Hubbell 3762C)	 Wall Receptacle <b>Hubbell CS6369</b> (Or Inline Hubbell CS6365C)

Install the Hubbell receptacle parallel to the wall  
and at least 24 inches from the floor.

- Power Consumption - Entire System 10 kVA (ISE, c503, and e801 plus sample supply unit and sample buffers)
- Power outlet (208V 50 amp) must be within 30 ft of the analyzer.



☐ Power outlet (115V 15 amp) must be within 8 ft of the printer/firewall/cobas link cart.

☐ It is recommended for the power to be fed from Emergency power (Generator backed up).

☐ A ground lead measuring less than 10 ohm with respect to earth ground is required. Safety ground must not be able to assume any potential above any grounded fixture located in the immediate vicinity of the analyzer (0.1 VAC maximum).

**NOTE:** Roche through an agreement with a power supply vendor provides an Uninterruptible Power Supply or Supplies (UPS) as part of the **cobas pro** analyzer series shipment. These UPS(s) provide power conditioning and short term AC power (5-10 minutes). If the site already has an Uninterruptible Supply built-into the facility, then it is not recommended that the Roche provided supplies be utilized (Uninterruptible supply should not be powered by another uninterruptible supply).

## ISE module

## c 503

## E 801

## Liquid Waste Volume

<input type="checkbox"/>	Concentrated liquid waste volume: < 2.0 L/h	Concentrated liquid waste volume: < 2 L/h	Concentrated liquid waste volume: ≤ 7 L/h
		Normal liquid waste volume: < 30 L/h	Normal liquid waste volume: < 22 L/h

☐ Floor drain required ≥ 2 inches ID and within 4 inches of the floor

☐ Floor drain must be located within 16 ft of the analyzer

☐ Waste is by gravity discharge

## Deionized Water Consumption

☐ Average in Operation: 2 L/h 32 L/h 30 L/h

☐ Clinical Laboratory Reagent Water (CLRW)

- Bacteria <10 CFU/mL Resistivity >10 MΩ.cm Total Organic Carbon(TOC) 500 PPB Particles 0.2um filtration or better

☐ Pressure 7.3 to 49.3 psi. Customer supplies a gauge and regulator valve to monitor and adjust water pressure.

☐ Temperature 54 to 86 °F (10 – 30°C)

☐ Customer supplies an approved shutoff valve near the analyzer. 5/8 hose barb with an internal diameter of 1/2 inch supplied by customer. DI water components must be either PVC or Stainless steel.

☐ Any dead leg tubing **must** be < 10 feet.

**Note:** It has been found that systems utilizing a re-circulation/polishing loop, with a minimal velocity of 5 ft/sec. and a 0.2 micron output filter, have been extremely effective in minimizing bacterial growth and maintaining overall water quality.

**Note:** If an acceptable floor or wall drain cannot be provided, the facility may provide, at their own discretion, an automated waste pumping solution (ie. separate waste container and "sump pump") capable of safely and effectively eliminating analyzer waste at the specified discharge rate.

## SSU

## ISE module

## c 503

## e 801

## Environmental

☐ Ambient Temperature 64 - 89° F in operation ≤ 3.6° F Δ/h

☐ Relative Humidity 30 – 85% non-condensing

☐ Noise emissions to Environment <65 dbA in Operation/<56 dbA in STBY

☐ Heat Dissipation

(BTU/h)w/SB

2730

1365

8190

6824



- |   |         |         |
|---|---------|---------|
| <input type="checkbox"/> Normal altitude above sea level  | 6562 ft | 6562 ft |
| <input type="checkbox"/> Max. altitude above sea level with high altitude kit   | 9843 ft | 9843 ft |
| <input type="checkbox"/> The environment should be relatively dust and vibration free.  |         |         |
| <input type="checkbox"/> The cobas pro integrated solutions should be located ample distance away from any equipment generating electromagnetic noise or electromagnetic wave interference (such as centrifuge, electric discharge machine, mobile telephone, transceiver, cordless telephone, etc) |         |         |

### Routing Considerations (Uncrated)

- ☐ Is there a loading dock or other suitable facility to allow the analyzer to be safely unloaded?
- ☐ Truck with lift gate or Forklift needed? Is a forklift available onsite?
- ☐ Is there any restriction to the length of truck that will be used to deliver the analyzer?
- ☐ Is there any restriction on delivery times/days?
- ☐ Is there an unpacking location available?

### Additional Considerations

- ☐ Installation Clearance: Minimum door width (opening) – 36 in. Minimum door height – 69 in. Minimum turning radius – 60 in.
- ☐ Elevator Clearance (if necessary): Width - 36 in. Depth - 48 in.
- ☐ Clearance for proper circulation and accessibility: Right: 24 in. Front: 39 in. Back: 28 in.
- ☐ Adequate ventilation in the immediate area of operation should be provided without airflow directly onto or across the top of the analyzer.
- ☐ The **cobas** pro analyzer should be located ample distance away from any equipment generating electromagnetic noise or electromagnetic wave interference such as centrifuge, electric discharge machine.
- ☐ System should be protected from direct sunlight.
- ☐ Adequate refrigerated storage at 2 to 8 °C **must** be available.
- ☐ Access for Maintenance is as follows:  
It is recommended to have **at least** 16" on the right side, 40" at the front and a necessary space of 28" at the back to ensure proper air circulation and accessibility for maintenance and operation.

### Network Connections

- ☐ A wired 10/100/1000 RJ45 female connection (1 drop per system) to the customer's network
- ☐ Roche provides a Firewall. The use of the Firewall is mandatory and requires a static IP address. See separate Firewall pre-site document.
- ☐ The communication is outbound through port 443 (HTTPS) and port 80 (HTTP with secure payload). Secure payload means that the data is encrypted, signed, and zipped. Protocols used are HTTP, HTTPS, SSH, SOAP, and BITS.
- ☐ For outbound LIS communication the Firewall will require access to the customers LIS over specific ports.
- ☐ For inbound LIS communication, the LIS will require access to the Firewall's static IP address.
- ☐ If installing at a VA account please follow IPB for Axeda setup – Veterans Affairs VPN Procedure
- ☐ For Roche remote diagnostics, Axeda, the Firewall's static IP address will need access to:

IP Address	Outbound Port	Host name
196.3.50.39	80	teleservice.roche.com**
62.209.44.11	443	remoteservice.roche.com
209.202.167.21	443	remoteservice-dr.roche.com
62.209.44.21	443	remoteservice-gas1.roche.com
62.209.44.22	443	remoteservice-gas2.roche.com
209.202.167.19	443	remoteservice-gas3.roche.com
209.202.167.20	443	remoteservice-gas4.roche.com
120.136.45.231	443	remoteservice-gas5.roche.com
120.136.45.230	443	remoteservice-gas6.roche.com



## Cart

- ☐ A cart will be **optional**. There are two options: 24 in wide standing cart and a 36 in wide sitting cart.
- ☐ Cart **must** be on the left side of the instrument under the monitor

Configurations	Total length Including Monitor (To stat port)
	in
<ISE   c503>	129 (113)
<ISE   c503   e801>	188 (172)
<ISE   c503   ISE   c503>	206 (190)
<ISE   c503   ISE   c503   e801>	265 (249)
<ISE   c503   e801   e801>	247 (231)
<ISE   c503   ISE   c503   e801   e801>	324 (308)
<ISE   c503   e801   e801   e801>	307 (291)
<e801>	111 (95)
<e801   e801>	171 (155)
<e801   e801   e801>	230 (214)
<e801   e801   e801   e801>	289 (273)




## Visual Guide ISE / c 503 / e 801



cobas pro  
19,109 BTUs  
UPS  
4,239 BTUs  
cobas Link  
132 BTUs  
Firewall  
51 BTUs  
Printer  
1877 BTUs

**Total BTUs  
25,408**





The cobas pro is approximately  
 (Not including UPS or computer table)  
**4,500 lbs**

# Visual Guide ISE / c 503 / e 801

## Water

CLRW grade DI water @ 7.3 to 50 PSI

Operational Consumption = 61 L/Hr  
Maximum Consumption = 64 L/Hr  
Required Flow Rate = 70 L/Hr






# Visual Guide

## ISE / c 503 / c 503 / e 801 / e 801



cobas pro  
35,488 BTUs  
2 UPS  
8,478 BTUs  
cobas Link  
132 BTUs  
Firewall  
51 BTUs  
Printer  
1877 BTUs

**Total BTUs**  
**46,026**



The cobas pro is approximately  
(Not including UPSs or computer table)  
**8,300 lbs**







# Visual Guide

## ISE / c 503 / c 503 / e 801 / e 801

### Water

CLRW grade DI water @ 7.3 to 50 PSI

Operational Consumption = 124 L/Hr  
 Maximum Consumption = 132 L/Hr  
 Required Flow Rate = 160 L/Hr



**Important:** The strict adherence to the site requirements contained in this document is an essential component in achieving a successful and compliant installation. Any delays in site requirement implementation over the jointly agreed dates can have a negative impact on and delay instrument shipment and/or installation. Deviations from these site specifications should be reported to the Roche Technical Support Center (1-800-428-2336).



This image shows a full page of blank graph paper. The grid consists of small, equal-sized squares formed by thin black lines. There are no margins, text, or other markings on the page.